

Made for
iPod iPhone iPad

Price List w.e.f. 1st March 2018

Introducing
Enya IIC



ReSound GN

ReSound Assist

More

hearing care
wherever you are.

ReSound Batteries



For better performance use
ReSound Batteries

Buy any 2 Wireless Hearing Aids of MRP 64,995/- each or higher,
Get any 1 Wireless Accessory worth MRP 17,700/-

FREE

EMI

Options Available

GN Making Life Sound Better

Smart Hearing 

All hearing aids are iSolate™ nanotech treated with 2 year warranty

Sr. No.	Model Name	No of Channel	Full on Gain 50dB	Battery Type	Fitting Range	Category	MRP Per Unit (Rs)
---------	------------	---------------	-------------------	--------------	---------------	----------	-------------------

PEDIATRIC PRODUCTS

Up Smart 9 (Made for iPhone)

1	Up Smart 998 DLW SP BTE	17 Channel	86	675	60-120	Profound	2,19,995/-
2	Up Smart 988 DLW P/HP BTE (Metal Hook)	17 Channel	74/80	13	25-105/115	Mild-Profound	

Up Smart 7 (Made for iPhone)

3	Up Smart 798 DLW SP BTE	17 Channel	86	675	60-120	Profound	1,09,995/-
4	Up Smart 788 DLW P/HP BTE (Metal Hook)	17 Channel	74/80	13	25-105/115	Mild-Profound	

Up Smart 5 (Made for iPhone)

5	Up Smart 598M DLW SP BTE	9 Channel	86	675	60-120	Profound	64,995/-
6	Up Smart 588M DLW P/HP BTE (Metal Hook)	9 Channel	74/80	13	25-105/115	Mild-Profound	

SUPER POWER PRODUCTS

ENZO 3D 9 (Made for iPhone + Remote Assist)

7	ENZO 3D 998 DW SP BTE	17 Channel	86	675	60-120	Profound	2,30,995/-
8	ENZO 3D 988 DW P/HP BTE (Metal Hook)	17 Channel	74/80	13	30-105/115	Mild-Profound	

ENZO 3D 7 (Made for iPhone + Remote Assist)

9	ENZO 3D 798 DW SP BTE	14 Channel	86	675	60-120	Profound	1,15,495/-
10	ENZO 3D 788 DW P/HP BTE (Metal Hook)	14 Channel	74/80	13	30-105/115	Mild-Profound	

ENZO 3D 5 (Made for iPhone + Remote Assist)

11	ENZO 3D 598 DW SP BTE	12 Channel	86	675	60-120	Profound	67,995/-
12	ENZO 3D 588 DLW P/HP BTE (Metal Hook)	12 Channel	74/80	13	30-105/115	Mild-Profound	

ENZO 5

13	ENZO 598 DW SP BTE (Non Mfi)	9 Channel	86	675	60-120	Profound	37,995/-
----	------------------------------	-----------	----	-----	--------	----------	----------

Magna

14	Magna 490 DVI SP BTE	9 Channel	86	675	60-120	Profound	28,495/-
15	Magna 290 DVI SP BTE	6 Channel	86	675	60-120	Profound	19,995/-

ReSound LiNX 3D™

LiNX 3D RECHARGEABLE (Rechargeable kit to be purchased separately)

16	LiNX3D 961 DRWZ RIE	17 Channel	61/67/74/82	312	20-95/110	Mild-Profound	2,62,495/-
17	LiNX 3D 761 DRWZ RIE	14 Channel	61/67/74/82	312	20-95/110	Mild-Profound	1,49,995/-
18	LiNX 3D 561 DRWZ RIE	12 Channel	61/67/74/82	312	20-95/110	Mild-Profound	86,495/-
19	Complete Rechargeable Kit (2 Rechargeable Batteries, 2 Rechargeable Battery Doors, 1 Recharging Dock)						23,600/-



LiNX 3D
Rechargeable Kit



LiNX 3D
Rechargeable

NOTE: Products (Rows) highlighted in blue colour are wireless products.

All prices are subject to change without prior notice.

All hearing aids are iSolate™ nanotech treated with 2 year warranty

Sr. No.	Model Name	No of Channel	Full on Gain 50dB	Battery Type	Fitting Range	Category	MRP Per Unit (Rs)
LiNX 3D 9 (Made for iPhone + ReSound Assist)							
20	LiNX3D 961 DW RIE	17 Channel	61/67/74/82	312	20-95/110	Mild-Profound	2,62,495/-
21	LiNX3D 962 DW RIE	17 Channel	61/67/74/82	13	20-95/110	Mild-Profound	
22	LiNX3D 988 DW P BTE (Metal Hook)	17 Channel	74	13	30-105	Mild-Profound	
23	LiNX3D 977 DW BTE Open/Classic	17 Channel	62/66	13	20-85/110	Mild-Profound	
24	LiNX3D 930 DW ITC/DW P ITC	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
25	LiNX3D 930 DW ITC/DW P ITC (Wireless CIC)	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	2,52,495/-
26	LiNX3D 930 D ITC/DP ITC	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
27	LiNX 3D 910 CIC/P CIC	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	
28	LiNX 3D 9 IIC	17 Channel	49(LP)	10A	20-85	Mild-Severe	2,62,495/-
LiNX 3D 7 (Made for iPhone + ReSound Assist)							
29	LiNX 3D 761 DW RIE	14 Channel	61/67/74/82	312	20-95/110	Mild-Profound	1,49,995/-
30	LiNX 3D 762 DW RIE	14 Channel	61/67/74/82	13	20-95/110	Mild-Profound	
31	LiNX 3D 788 DW P BTE (Metal Hook)	14 Channel	74	13	30-105	Mild-Profound	
32	LiNX 3D 777 DW BTE Open / Classic	14 Channel	62/66	13	20-85/100	Mild-Profound	
33	LiNX 3D 730 DW ITC/DW P ITC	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
34	LiNX 3D 730 W ITC/W P ITC (Wireless CIC)	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
35	LiNX 3D 730 D ITC/DP ITC	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	1,39,995/
36	LiNX 3D 710 CIC/P CIC	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	1,49,995/-
37	LiNX 3D 7 IIC	14 Channel	49(LP)	10A	20-85	Mild-Severe	
LiNX 3D 5 (Made for iPhone + ReSound Assist)							
38	LiNX 3D 561 DW RIE	12 Channel	61/67/74/82	312	20-95/110	Mild-Profound	86,495 /-
39	LiNX 3D 562 DW RIE	12 Channel	61/67/74/82	13	20-95/110	Mild-Profound	
40	LiNX 3D 588 DW P BTE (Metal Hook)	12 Channel	74	13	30-105	Mild-Profound	
41	LiNX 3D 577 DW BTE Open / Classic	12 Channel	62/66	13	20-85/100	Mild-Profound	
42	LiNX 3D 530 DW ITC/DW P ITC	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
43	LiNX 3D 530 W ITC/W P ITC (Wireless CIC)	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	78,995 /-
44	LiNX 3D 530 D ITC/DP ITC	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
45	LiNX 3D 510 CIC/P CIC	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	
46	LiNX 3D 5 IIC	12 Channel	49(LP)	10A	20-85	Mild-Severe	86,495/-

NOTE: Products (Rows) highlighted in blue colour are wireless products.


All prices are subject to change without prior notice.

All hearing aids are iSolate™ nanotech treated with 2 year warranty


Sr. No.	Model Name	No of Channel	Full on Gain 50dB	Battery Type	Fitting Range	Category	MRP Per Unit (Rs)
---------	------------	---------------	-------------------	--------------	---------------	----------	-------------------

ReSound Enya™


ENYA 4

46	Enya 462 DW RIE	10 Channel	61/67/74/82	312	20-90/110	Mild-Profound	49,995/-
47	Enya 488 DW P/HP BTE (Metal Hook)	10 Channel	74/80	13	25-105/115	Mild-Profound	
48	Enya 477 DW BTE Open / Classic	10 Channel	62/66	13	10-80/100	Mild-Profound	
49	Enya 430 DW ITC/DW P ITC	10 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
50	Enya 430 D ITC/DP ITC	10 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	49,995/-
51	Enya 410 CIC/P CIC	10 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	
52	Enya 410 CIC (IIC with black face plate) 	10 Channel	49(LP)	10A	20-85	Mild-Severe	

ENYA 3

53	Enya 362 DW RIE	8 Channel	61/67/74/82	312	20-90/110	Mild-Profound	34,995/-
54	Enya 388 DW P/HP BTE (Metal Hook)	8 Channel	74/80	13	25-105/115	Mild-Profound	
55	Enya 377 DW BTE Open / Classic	8 Channel	62/66	13	10-80/100	Mild-Profound	
56	Enya 330 DW ITC/DW P ITC	8 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
57	Enya 330 D ITC/DP ITC	8 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	34,995/-
58	Enya 310 CIC/P CIC	8 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	
59	Enya 310 CIC (IIC with black face plate) 	8 Channel	49(LP)	10A	20-85	Mild-Severe	

ENYA 2

60	Enya 288 DW P/HP BTE (Metal Hook)	6 Channel	74/80	13	25-105/115	Mild-Severe	24,995/-
61	Enya 277 DW BTE Open / Classic	6 Channel	62/66	13	10-80/100	Mild-Profound	
62	Enya 230 DW ITC/DW P ITC	6 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
63	Enya 230 D ITC/DP ITC	6 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	24,995/-
64	Enya 210 CIC/PCIC	6 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	
65	Enya 210 CIC (IIC with black face plate) 	6 Channel	49(LP)	10A	20-85	Mild-Severe	

Sr. No.	Model Name	Features	Full on Gain 50dB	Battery Type	Fitting Range	MRP perUnit (Rs)
---------	------------	----------	-------------------	--------------	---------------	------------------

ReSound Vea™

Vea 3

66	Vea 362 DVI RIE	<ul style="list-style-type: none"> • 7 Band & 7 Channel WARPTM • 2 Environmental Programs + Telecoil + DAL • Adaptive Directionality™, SoftSwitching™ • Fixed Directionality • Noise Tracker™ II • Dual Stabilizer® II DFS with Whistle Control™ • Auto Phone Detection: Phone Now • Data Logging • In-situ Audiometry 	61/69	312	10-85/100	21,995 /-
67	Vea 388 DVI POWER BTE		75	13	25-115	29,995 /-
68	Vea 377 DVI BTE Open / Classic		69	13	10-80/110	28,995 /-
69	Vea 330 D ITC/DP ITC		52/60	312	15-80/20-90	29,995/-
70	Vea 310 CIC/P CIC		44/52	10A	10-80/90	

NOTE: Products (Rows) highlighted in blue colour are wireless products.

All prices are subject to change without prior notice.

All hearing aids are iSolate™ nanotech treated with 2 year warranty

Sr. No.	Model Name	Features	Full on Gain 50dB	Battery Type	Fitting Range	MRP perUnit (Rs)
Vea 2						
72	Vea 280 DVI Power BTE	<ul style="list-style-type: none"> • 6 Band & 6 Channel WARP™ • 2 Environmental Programs + Telecoil/DAI • Fixed Directionality • Noise Tracker™ II • Dual Stabilizer® II DFS with WhistleControl™ • Auto Phone Detection : Phone Now • Data Logging • In-situ Audiometry 	75	13	25-115	18,995 /-
73	Vea 270 DVI BTE Open / Classic		69	13	10-80/110	17,995 /-
74	Vea 230 D ITC/DP ITC		52/60	312	15-80/20-90	19,995 /-
75	Vea 210 D ITC/DP ITC		44/52	10A	10-18/90	

Vea 1						
76	Vea 180 VI Power BTE	<ul style="list-style-type: none"> • 4 Band & 4 Channel WARP™ • 1 Environment Programs + Telecoil + DAI • Noise Tracker™ • Stabilizer® II DFS with Whistle Control™ • Data Logging • In-situ Audiometry 	75	13	25-115	13,995 /-
77	Vea 170 VI BTE Open / Classic		69	13	10-80/110	12,995 /-
78	Vea 130 D ITC/DP ITC		52/60	312	15-80/20-90	14,995 /-
79	Vea 110 ITC/DP ITC		44/52	10A	10-18/90	

Sr. No.	Model Name	No of Channel	Full on Gain 50dB	Battery Type	Fitting Range	Category	MRP Per Unit (Rs)
---------	------------	---------------	-------------------	--------------	---------------	----------	-------------------

ReSound ENZO²

ENZO² 9 (Made for iPhone)

80	ENZO ² 998 DW SP BTE	17 Channel	86	675	60 - 120	Profound	2,19,995 /-
81	ENZO ² 988 DLW P/HP BTE (Metal Hook)	17 Channel	74/80	13	30-105/115	Mild-Profound	

ENZO² 7 (Made for iPhone)

82	ENZO ² 798 DW SP BTE	14 Channel	86	675	60-120	Profound	1,09,995 /-
83	ENZO ² 788 DLW P/HP BTE (Metal Hook)	14 Channel	74/80	13	30-105/115	Mild-Profound	

ENZO² 5 (Made for iPhone)

84	ENZO ² 598 DW SP BTE	12 Channel	86	675	60 - 120	Profound	64,995 /-
85	ENZO ² 588 DLW P/HP BTE (Metal Hook)	12 Channel	74/80	13	30-105/115	Mild-Profound	

ReSound LiNX²™

LiNX² 9 (Made for iPhone)

86	LiNX ² 961 DW RIE	17 Channel	61/67/74/82	312	20-95/110	Mild-Profound	2,49,995/-
87	LiNX ² 962 DW RIE	17 Channel	61/67/74/82	13	20-95/110	Mild-Profound	
88	LiNX ² 988 DW P BTE (Metal Hook)	17 Channel	74	13	30-105	Mild-Profound	
89	LiNX ² 977 DW BTE Open / Classic	17 Channel	62/66	13	20-85/100	Mild-Profound	
90	LiNX ² 930 DW ITC/DW P ITC	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/100	Mild-Profound	
91	LiNX ² 930 W ITC/W P ITC (Wireless CIC)	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/100	Mild-Profound	
92	LiNX ² 930 D ITC/DP ITC	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/100	Mild-Profound	
93	LiNX ² 910 CIC/P CIC	17 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/100	Mild-Profound	
94	LiNX ² 9 IIC	17 Channel	49(LP)	10A	20-85	Mild-Severe	

NOTE: Products (Rows) highlighted in blue colour are wireless products.

All prices are subject to change without prior notice.

All hearing aids are iSolate™ nanotech treated with 2 year warranty

Sr. No.	Model Name	No of Channel	Full on Gain 50dB	Battery Type	Fitting Range	Category	MRP Per Unit (Rs)
LiNX² 7 (Made for iPhone)							
95	LiNX ² 761 DW RIE	14 Channel	61/67/74/82	312	20-95/110	Mild-Profound	1,29,995/-
96	LiNX ² 762 DW RIE	14 Channel	61/67/74/82	13	20-95/110	Mild-Profound	
97	LiNX ² 788 DW P BTE (Metal Hook)	14 Channel	74	13	30-105	Mild-Profound	
98	LiNX ² 777 DW BTE Open / Classic	14 Channel	62/66	13	20-85/100	Mild-Profound	
99	LiNX ² 730 DW ITC/DW P ITC	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
100	LiNX ² 730 W ITC/W P ITC (Wireless CIC)	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	1,14,995/-
101	LiNX ² 730 D ITC/DP ITC	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
102	LiNX ² 710 CIC/P CIC	14 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	1,29,995/-
103	LiNX ² 7IIC	14 Channel	49(LP)	10A	20-85	Mild-Severe	

LiNX² 5 (Made for iPhone)							
104	LiNX ² 561 DW RIE	12 Channel	61/67/74/82	312	20-95/110	Mild-Profound	82,995/-
105	LiNX ² 562 DW RIE	12 Channel	61/67/74/82	13	20-95/110	Mild-Profound	
106	LiNX ² 588 DW P BTE (Metal Hook)	12 Channel	74	13	30-105	Mild-Profound	
107	LiNX ² 577 DW BTE Open / Classic	12 Channel	62/66	13	20-85/100	Mild-Profound	
108	LiNX ² 530 DW ITC/ DW P ITC	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
109	LiNX ² 530 W ITC/ W P ITC (Wireless CIC)	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	74,995/-
110	LiNX ² 530 D ITC/ DP ITC	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	312	20-85/110	Mild-Profound	
111	LiNX ² 510 CIC/P CIC	12 Channel	49(LP)/59 (MP) 69(HP)/79 (UP)	10A	20-85/110	Mild-Profound	82,995/-
112	LiNX ² 5 IIC	12 Channel	49(LP)	10A	20-85	Mild-Severe	

Sr. No.	Model Name	Features	Full on Gain 50dB	Battery Type	Fitting Range	MRP perUnit (Rs)
---------	------------	----------	-------------------	--------------	---------------	------------------

Trimmer Digitals

Match 3 (2 Channels)

113	Match 3T 90 VI BTE	<ul style="list-style-type: none"> • 2 Band & 2 Channel WARP™ Trimmer Digital • 4 Programs (2 Environmental Programs +Telecoil + DAI) • Noise Reduction • Feedback manager • L-Trimmer, H-Trimmer, P-Trimmer 	86	675	60-120	13,095 /-
114	Match 3T 80 V POWER BTE		72	13	25-110	9,995 /-
115	Match 3T 70 V BTE		64	13	20-90	9,445 /-

Match 2 (2 Channels)

116	Match 2T 80 V POWER BTE	<ul style="list-style-type: none"> • 2 Band & 2 Channel WARP™ Trimmer Digital • 2 Environmental Programs <ul style="list-style-type: none"> • Noise Reduction • Feedback Manager • L-Trimmer, P-Trimmer 	72	13	25-110	8,895 /-
117	Match 2T 70 V BTE		64	13	20-90	8,395 /-
118	Match 2T 50 P ITE		68	13	20-95	11,995 /-
119	Match 2T 30 P ITC		60	312	20-85	
120	Match 2T 10 P CIC		52	10A	20-75	

All prices are subject to change without prior notice.

All hearing aids are iSolate™ nanotech treated with 2 year warranty

Sr. No.	Model Name	Features	Full on Gain 50dB	Battery Type	Fitting Range	MRP per Unit (Rs)
Match 1 (2 Channels)						
121	Match 1T 70 V BTE	<ul style="list-style-type: none"> • 2 Band & 2 Channel WARP™ Trimmer Digital • 2 Environmental Programs • Noise Reduction • Feedback Manager • L-Trimmer 	64	13	20-90	7,495 /-
122	Match 1T 50 P ITE		68	13	20-95	
123	Match 1T 30 P ITC		60	312	20-85	9,995 /-
124	Match 1T 10 P CIC		52	10A	20-75	
Analog BTE						
125	215 UT	<ul style="list-style-type: none"> • Analog • M-T-O Switch • H-Trimmer • Telecoil 	66	675	30-90	7,495 /-
126	793 U HG	<ul style="list-style-type: none"> • Analog • N-H-O Switch 	52	13	20-70	5,745 /-

Accessories MRP Inclusive of GST

Sr.No.	LINX 3D - Rechargeable Accessories		MRP per Unit (Rs)
1	Rechargeable Battery Door		6,077/-
2	Rechargeable Battery		3,136/-
3	Rechargeable Dock		6,313/-
Sr.No.	Up Smart Accessories		MRP per Unit (Rs)
4	BEND and GO		1,000/-
5	STICK and GO		1,000/-
6	CLIP and GO		600/-
7	Pediatric Care Kit		6,000/-
Sr.No.	Unite Wireless Accessories (with 1 year warranty)		MRP per Unit (Rs)
8	PHONE CLIP +		17,700/-
9	TV UNITE 2		
10	REMOTE CONTROL 2		
11	MULTI MIC		
Sr.No.	Accessories		MRP Per
12	RIE RECEIVERS	AVAILABLE IN 0,1, 2 & 3 SIZES	10,000/-
13	MICROMOULD FOR RIE PRODUCTS		600/-
14	HARD BTE EARMOULD		600/-
15	SOFT BTE EARMOULD		900/-
16	HF3 FILTER	STRIP INCLUDES 15 FILTER PIECES	1,050/-
17	CERUSTOP WAXGUARD BOX	STRIP INCLUDES 8 FILTER PIECES	700/-
18	RESHELLING CHARGES	CHANGE IN CANAL SIZE, DAMAGED ETC.	4,000/-
Sr.No.	Extended Warranty		MRP Per
19	1 YEAR EXTENDED WARRANTY	APPLICABLE ON MRP ≥ RS. 15,000/-	10% of MRP + 18% GST
20	2 YEAR EXTENDED WARRANTY		15% of MRP + 18% GST



MULTI MIC



TV UNITE 2



PHONE CLIP +



REMOTE CONTROL 2

NOTE: Products (Rows) highlighted in blue colour are wireless products.

All prices are subject to change without prior notice.

FEATURES

World's Smartest Hearing Aids

BENEFITS

What more you get with LiNX 3D 9

- + Binaural Directionality III
- + Spatial Sense
- + Sound Enhancer in Smart 3D App

LiNX 3D 9

17 Channel

- + With BD-III you can focus on the sounds that are important to you without losing touch with what's going on around you. Providing the brain with the right information from both ears lets you effortlessly follow conversations and shift attention naturally.
- + Supports your ability to detect where sounds are coming from, so you experience a more natural sense of your surroundings and the best possible sound quality.
- + Provides you full control of what you want to hear, where you want to hear from. Using your smartphone you can make your hearing aid focus on the person you want to hear in a group of people & control the noise from the surroundings

What more you get with LiNX 3D 7

- + Binaural Directionality
- + Environmental Optimizer

LiNX 3D 7

14 Channel

- + Both your hearing aids automatically focus on the sounds you're most interested in while maintaining a sense of what is going on around you. This supports your natural ability to focus on speech, even in the noisiest of situations.
- + You can enjoy sound quality and listening comfort even when you are rapidly changing sound environments, like from office to outdoor or party to traffic etc, automatically without a need to press a button

✓ ReSound Assist

✓ Smart Range II Chip

✓ Noise Tracker II

✓ Made for All

✓ Smart 3D App

✓ Connects to Unite Wireless Accessories

LiNX 3D 5

12 Channel

- ✓ Sometimes the hearing aid settings you like in the clinic perform differently in your everyday life. If you're experiencing recurring problems and need more fine-tuning sessions you need not go to hearing clinic, your hearing aids can be programmed wherever you are. "Yes you can get your hearing aids fine-tuned by your hearing care professional, without going to hearing clinic", using Resound Assist in Smart 3D App.
- ✓ Faster sound processing, decision making by the processor chip makes hearing feel effortless. Once you are home after a full day, you'll still have energy to spare for your loved ones. Because all that hearing, which used to drain your energies, now seems so effortless, so natural.
- ✓ Improves your comfort in noisy listening situations like restaurants, railway stations without compromising speech understanding.
- ✓ Your hearing aids function like bluetooth headphones for making phone calls, streaming music, watching videos etc.
- ✓ ReSound LiNX 3D automatically adjusts to your sound environment. Still, in case you want to personalize and control your sound at any time, you can do it with the new ReSound Smart 3D app and save it as a Favorite program. You can even geo-tag it to automatically activate when you arrive at work, home, the golf course etc.
- ✓ Accessory (TV Streamer 2) plugged into your TV, can stream stereo sound directly into your hearing aids, at your preferred volume while the people around you listen at their comfortable level with out disturbing each other.

FEATURES

Smart Super Power Hearing Aids

BENEFITS

What more you get with ENZO 3D 9

- + Binaural Directionality III
- + Spatial Sense
- + Sound Enhancer in Smart 3D App

ENZO 3D 9

17 Channel

- + With BD-III you can focus on the sounds that are important to you without losing touch with what's going on around you. Providing the brain with the right information from both ears lets you effortlessly follow conversations and shift attention naturally.
- + Supports your ability to detect where sounds are coming from, so you experience a more natural sense of your surroundings and the best possible sound quality.
- + Provides you full control of what you want to hear, where you want to hear from. Using your smartphone you can make your hearing aid focus on the person you want to hear in a group of people & control the noise from the surroundings

What more you get with ENZO 3D 7

- + Binaural Directionality
- + Environmental Optimizer

ENZO 3D 7

14 Channel

- + Both your hearing aids automatically focus on the sounds you're most interested in while maintaining a sense of what is going on around you. This supports your natural ability to focus on speech, even in the noisiest of situations.
- + You can enjoy sound quality and listening comfort even when you are rapidly changing sound environments, like from office to outdoor or party to traffic etc, automatically without a need to press a button

✓ ReSound Assist

✓ Smart Range II Chip

✓ DFS-Ultra II With Music mode

✓ Noise Tracker II

✓ Made for All

✓ Smart 3D App

✓ Connects to Unite Wireless Accessories

ENZO 3D 5

12 Channel

- ✓ Sometimes the hearing aid settings you like in the clinic perform differently in your everyday life. If you're experiencing recurring problems and need more fine-tuning sessions you need not go to hearing clinic, your hearing aids can be programmed wherever you are. "Yes you can get your hearing aids fine-tuned by your hearing care professional, without going to hearing clinic", using Resound Assist in Smart 3D App.
- ✓ Faster sound processing, decision making by the processor chip makes hearing feel effortless. Once you are home after a full day, you'll still have energy to spare for your loved ones. Because all that hearing, which used to drain your energies, now seems so effortless, so natural.
- ✓ These are powerful hearing aids and there is a need to control this power effectively, otherwise it leads to whistling sound, DFS enabled with adaptive filter never lets your hearing aids whistle. Music Mode added to DFS ultra-II differentiates musical tones from whistle and cancels only whistling sound from the hearing aids.
- ✓ Improves your comfort in noisy listening situations like restaurants, railway stations without compromising speech understanding.
- ✓ Your hearing aids function like bluetooth headphones for making phone calls, streaming music, watching videos etc.
- ✓ ReSound ENZO 3D automatically adjusts to your sound environment. Still, in case you want to personalize and control your sound at any time, you can do it with the new ReSound Smart 3D app and save it as a Favorite program. You can even geo-tag it to automatically activate when you arrive at work, home, the golf course etc.
- ✓ Accessory (TV Streamer 2) plugged into your TV, can stream stereo sound directly into your hearing aids, at your preferred volume while the people around you listen at their comfortable level with out disturbing each other.

FEATURES

Advanced Hearing Aids

BENEFITS

What more you get with Enya 4

- + Ear to Ear Communication
- + Comfort Phone

Enya 4
10 Channel

- + You can change the volume or program of both hearing aids by adjusting just one of the hearing aids, thanks to Ear-to-Ear communication.
- + When to receive a phone call in a noisy situation, your Enya 4 hearing aids communicate to each other to provide you best phone reception in one ear and the other hearing aid allowing you to be aware of environmental sounds.

What more you get with Enya 3

- + Soft Switching
- + Control App & Tinnitus Relief App
- + Wind Guard
- + Connects to Unite Wireless Accessories

Enya 3
8 Channel

- + With this feature hearing aids automatically choose to bring a balance between better understanding of speech in noise and awareness of the environment. i.e. when you move to noisy situation hearing aid helps you understand speech and also be aware of warning sounds like vehicle horns and alarm bells etc.
- + These free apps allow you to control your hearing aids volume, change the programs of the hearing aids, listen to a wide variety of tinnitus reduction music.
- + Reduces the discomforting wind noise experienced when you drive a two wheeler, car with open windows, live near a breeze area like seashore etc.
- + Accessory (e.g.:TV-Streamer-2) plugged into your TV, can stream stereo sound directly into your hearing aids, at your preferred volume. Phone Clip + routes your phone calls from your phone directly into the hearing aids, making your hearing aids function as Bluetooth headphones.

- ✓ Noise Tracker II
- ✓ DFS Ultra
- ✓ High Power & Discreet Design
- ✓ Adaptive Directionality
- ✓ Tinnitus Sound Generator

Enya 2
6 Channel

- ✓ Improves your comfort in noisy listening situations like restaurants, railway stations without compromising speech understanding.
- ✓ Feedback or whistling sound is an annoying sound made by hearing aids when you wear hearing aids and get close to walls, people etc. DFS ultra enabled with adaptive filter never lets your hearing aids whistle, even when you wear a scarf, helmet close to your hearing aids.
- ✓ You can now be fitted with hearing aids completely in the canal which are least visible and still get all the power that you need from your hearing aids unlike other hearing aids which compromise on visibility.
- ✓ Identify and understand speech clearly when in noisy environments with multiple speakers.
- ✓ Produces Tinnitus reducing sounds to relieve you of annoying, hallucinating ringing sensations in your ears.

FEATURES

Advanced SP Hearing Aids

BENEFITS

- + Smart Range Chip
- + DFS-Ultra II With Music mode
- + Connects to Unite WirelessAccessories
- + Ear to Ear Communication

-
- + Soft Switching
 - + Adaptive Directionality With Integrated Wind Noise Cancellation
 - + DFS Ultra
 - + Sound Shaper

-
- ✓ Noise Tracker II
 - ✓ Dual Stabilizer II DFS with Whistle Control
 - ✓ High Power
 - ✓ Cosmetic Design

Enzo 5

9 Channel

Magna 4

9 Channel































Magna 2

6 Channel

- + Faster sound processing, decision making by the processor chip makes hearing feel effortless. Once you are home after a full day, you'll still have energy to spare for your loved ones. Because all that hearing, which used to drain your energies, now seems so effortless, so natural.
- + These are powerful hearing aids and there is a need to control this power effectively, otherwise it leads to whistling sound, DFS enabled with adaptive filter never lets your hearing aids whistle. Music Mode added to DFS ultra-II differentiates musical tones from whistle and cancels whistling sound from hearing aids only.
- + Accessory (TV Streamer 2) plugged into your TV, can stream stereo sound directly into your hearing aids, at your preferred volume while the people around you listen at their comfortable level with out disturbing each other.
- + You can change the volume or program of both hearing aids by adjusting just one of the hearing aids, thanks to Ear-to-Ear communication.

-
- + With this feature hearing aids automatically choose to bring a balance between better understanding of speech in noise and awareness of the environment. i.e. when you move to noisy situation hearing aid helps you understand speech and also be aware of warning sounds like vehicle horns and alarm bells etc.
 - + This advanced directionality feature makes it easier to distinguish where speech and other important sounds are coming from improving speech understanding in noise. Furthermore, ReSound Magna reduces wind noise, enabling you to hear more comfortably outdoors.
 - + Feedback or whistling sound is an annoying sound made by hearing aids when you wear hearing aids and get close to walls, people etc. DFS ultra enabled with adaptive filter never lets your hearing aids whistle, even when you wear a scarf, helmet close to your hearing aids.
 - + This advanced algorithm makes sounds that are normally impossible to hear, audible again. This means that important speech cues are enhanced e.g. sounds such as children's voices and bird songs are clearly audible again.

-
- ✓ Improves your comfort in noisy listening situations like restaurants, railway stations without compromising speech understanding.
 - ✓ Feedback or whistling sound is an annoying sound made by hearing aids when they are providing high amplification. This can be controlled to a great extent.
 - ✓ Individuals with severe hearing loss need adequate amplification not only to hear sound but understand it. Resound Magna is one the most powerful BTEs in the market, providing more amplification in the situations where you need it the most.
 - ✓ Smallest size in the segment.

	9 Premium 	7 Premium 	5 Premium 	4 & 3 Essential 	2 & 1 Essential 
WARP BANDS	17	14 Up Smart: 17	12 Up Smart: 9	10 Magna: 9 ENYA: 8 Ve: 7	6 Ve: 4
SOUND SHAPER	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/> Magna	<input type="checkbox"/>
SOUND ENHANCER SMART 3D APP	ULTIMATE 	ADVANCED 	ADVANCED 		
PERSONALIZED NOISE REDUCTION	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
NOISE TRACKER II	ULTIMATE 	ADVANCED 	STANDARD 	BASIC 	 Ve 1: Noise Tracker
WINDGUARD	ULTIMATE 	ADVANCED 	STANDARD 	BASIC 	 INTEGRATED WIND NOISE MANAGEMENT Not in Ve 1
DIRECTIONALITY	6 Options including BIN-III with spatial sense in LiNX 3D & ENZO 3D and BIN-II in LiNX ² & ENZO ² only	4 Options including Binaural Directionality in LiNX 3D, ENZO 3D, LiNX ² & ENZO ² only	4 Options including Natural Directionality in LiNX 3D, ENZO 3D, LiNX ² & ENZO ² only	4 OPTIONS including Soft Switching	3 OPTIONS Ve 2 : 2 (Fixed & Omni) Ve 1 : 1 (Omni)
ENVIRONMENTAL OPTIMIZER	BINAURAL EO II 	EO 			
SYNCHRONIZED ACCEPTANCE MANAGER	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
DIGITAL FEEDBACK SUPPRESSION	DFS ULTRA II 	DFS ULTRA II 	DFS ULTRA II 	DFS ULTRA II  Ve: Dual Stabilizer with whistle control	DFS ULTRA II  Ve 1: Dual Stabilizer II with whistle control Ve 2 & Magna: Dual Stabilizer II with whistle control
REMOTE FINE TUNING	<input checked="" type="checkbox"/> LiNX 3D ENZO 3D	<input checked="" type="checkbox"/> LiNX 3D ENZO 3D	<input checked="" type="checkbox"/> LiNX 3D ENZO 3D	<input type="checkbox"/>	<input type="checkbox"/>
FULLY FLEXIBLE PROGRAMS	4	4	4	4 Ve : 2	3 ENYA Ve 1 : 1, Ve 2:2
LOW FREQUENCY BOOST	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
PHONE NOW	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/> Except Ve 1
SYNCHRONIZED VOLUME CONTROL / PUSH BUTTON DEVICE WITH BUTTONS	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/> Only in ENYA 4	<input type="checkbox"/>
NUMBER OF GAIN HANDLES	Max: 17	Max: 14	Max: 12	10 ENYA 3: 8 Magna & Ve: 6	6 Ve: 4
RESOUND UNITE MADE FOR IPHONE APPS	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> Smart 3D App in LiNX 3D & ENZO 3D	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> Smart 3D App in LiNX 3D & ENZO 3D	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> Smart 3D App in LiNX 3D & ENZO 3D	<input checked="" type="checkbox"/> Only ENYA <input type="checkbox"/> <input checked="" type="checkbox"/> Control App only in ENYA	<input checked="" type="checkbox"/> Remote Control in ENYA only <input type="checkbox"/>
TINNITUS RELIEF APP	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/> Only ENYA	<input type="checkbox"/>
INBUILT TINNITUS SOLUTION	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/> Only ENYA	<input checked="" type="checkbox"/> Only ENYA

ReSound GN

GN Hearing India Private Limited

Office No - 1201 & 1203, 12th Floor, V Times Square, Plot No - 3, Sector - 15, Palm Beach Road, CBD - Belapur, Navi Mumbai - 400614 (MS)

For any Query give Missed Call on +91 22 2752 4700

customercare.india@gnresound.com • www.resound.com/en-IN